

Now you can go anywhere with the Trav. Compact and lightweight, the Trav makes a great supplement to your primary wheelchair. The Trav handles easily in narrows spaces. It is perfect for use in restaurants, theaters, traveling by car, bus, train, or plane. Best of all, the Trav is built to the same high standards if quality and safety as all Nissin wheelchairs, giving you the freedom to go more places while also enjoying your peace of mind.

The Trav....

\$1,472 RETAIL

-Weighs less than 14lbs. -13.5(small) 15"(regular) wide. -Folds easily with the touch of one button. -Folds to meet airline requirements for carry-on bags. -Carry bag included for convenient storage and transportation.



ONE-TOUCH OPERATION , ANYTIME, ANYWHERE



Colours is the official distributor of Nissin Medical



Quality engineered for comfort, mobility and safety.

Nissin engineered the Trav with the help from researchers at a major rehabilitation center in Japan with the needs and desires of their clients in mind.

-The rear wheels are installed inside the chairs frame to reduce its width, increase maneuverability, and yet provide the comfort of a full size wheelchair. They have a step break for the convenience of caregivers.

-The thigh belt ensures balance seating.

-The flip-up armrest facilitates smoother transfers from the chair's seat.

The Trav allows you to enjoy more mobility while feeling safe and secure.



Size	Small	Regular
Seat Size	13.5" width x 13" depth	15" width x 13" depth
Overall Width (open)	33" height x 15.75" width x 25.6" depth	33" height x 17" width x 25.6" depth
Overall Width (closed)	21" height x 15.75" width x 7.5" depth	21" height x 17" width x 7.5" depth
Seat to Floor Height	19" front, 18" rear	19" front, 18" rear
Wheel Size	3" front, 5" rear	3" front, 5" rear
Frame	Aluminum	Aluminum
Frame Color	Silver	Silver
Upholstery Color	Black	Black

Colours is the official distributor of Nissin Medical

> Shipping Dimensions: 9W x 18.3D x 22.4H Shipping Weight: 16.5lbs

> > Return Policy: 60 days limit

- Product must be unused, clean, and not damaged.
- -Original packaging and materials intact.
- -Return shipping costs are the responsibility of the customer.

To place your order contact Colours customer service 951-808-9131 or email customerservice@colourswheelchair.com



Colours In Motion 2061 Cecilia Circle Corona, CA, 92881 TEL: 951-808-9131 Toll Free: 1-800-892-8998

